

# SPEAKMAN®

## SE-4300 & SE-4320

### GRAVITYFLO® PORTABLE EYEWASH

#### DESCRIPTION

The SE-4300 portable, gravity operated eyewash is the next best thing to a plumbed eyewash. The 20 gallon capacity provides a run time in excess of the ANSI Z358.1 required 15 minutes and offers the best performing spray. A simple pull of an elastomeric strap activates the unit and provides hands free operation. The integral drain allows flushing fluids to be contained avoiding slip hazards.

The SE-4300 is the industry's only portable eyewash that offers an optional side fitting drench hose. The optional SE-4920 drench hose attaches easily to the side discharge port and provides the ability for spot body flushing.

Includes one bottle of SE-4210 water preservative.

#### SPECIFICATIONS:

**Tank:** 20 gallon capacity. High density molded polyethylene

**Sprays:** Molded plastic

**Activator:** Elastomeric pull strap

**Mounting:** Shelf or wall mount (bracket included)

**Fill Cap:** Plastic fill cap with integral tamper prevention

**Side Outlet:** Integral 1/2" NHT fitting for attachment of SE-4920 drench hose.

**Drain:** Integral basin drain with 1/2" NHT fitting

**Weight:** (shipping weight): 18 lbs., **Filled:** 184 lbs.

**Compliance:** ANSI Z358.1

**Signs:** Integral Eyewash Sign

#### ADDITIONAL MODELS

- SE-4320 Gravity-operated eyewash including the side fitting drench hose
- SE-4350 Heated gravity fed eyewash provides protection in unheated spaces
- SE-4380 SE-4300 with cart & fluid collection stand

#### Available Accessories

**SE-4920** Side fitting drench hose with mounting hardware. Retrofittable for existing SE-4300 Gravityflo units.

**SE-4210** Water preservative

(Case of 4 - 8 oz. bottles).

Treats 5-21 gallons for up to 180 days.

**SE-4240** Fluid collection station

**SE-4360** Cart to transport SE-4300 series gravity-fed eyewashes

Unit exceeds the minimum 0.4 gpm for 15 minutes as stated by ANSI



Job Name _____
Location _____
Architect/Engineer _____
Sales Representative _____
Wholesaler _____
Date _____ Model # _____